

# Getein208

Hand-held Integrated System

Mini in One Way,

**Mighty** in Every Way



## Dual methodologies

Dry clinical chemistry and Immunofluorescence  
Supporting multiple test items  
Combo test, highly efficient



### Getein208

Hand-held Integrated System



## Small in size, easy to carry

Dimension: 161x88x24 mm  
Weight: 256 g



## Multi-scene application

Ambulance, Emergency Dept.,  
Clinical departments, ICU, etc.



## Available to multiple sample types

Whole blood, serum, plasma, peripheral blood,  
urine, etc.

- **Single Packaging**

The test cards are individually packaged.  
Convenient to use, no contamination to others

- **Lithium-ion Battery**

Battery life supports up to 300 tests.  
More than two months standby time

- **IoT Chip**

LIS and bluetooth printer connectivity

- **Removable Cover Plate for Cleaning**

Preventing cross-contamination



# Test Items

Items with \* are coming soon.

| Test Items                    | Diseases                                          | Sample Types           |
|-------------------------------|---------------------------------------------------|------------------------|
| <b>Cardiac Markers</b>        |                                                   |                        |
| *TnT                          | Myocardial Infarction                             | S/P/WB                 |
| cTnI                          | Myocardial Infarction                             | S/P/WB                 |
| D-Dimer                       | Venous Thromboembolism                            | P/WB                   |
| BNP                           | Heart Failure                                     | P/WB                   |
| NT-proBNP                     | Heart Failure                                     | S/P/WB                 |
| NT-proBNP/cTnI                | Heart Failure                                     | S/P/WB                 |
| CK-MB/cTnI/Myo                | Myocardial Damage/Infarction                      | S/P/WB                 |
| hs-CRP                        | Cardiovascular Inflammation                       | S/P                    |
| <b>Inflammation</b>           |                                                   |                        |
| *IL-6                         | Acute Inflammation                                | S/P/WB/Fingertip Blood |
| PCT                           | Sepsis, Bacterial Infection                       | S/P/WB/Fingertip Blood |
| hs-CPR+CRP                    | Cardiovascular Inflammation, Normal Inflammation  | S/P/WB                 |
| PCT/CRP                       | Sepsis, Bacterial Infection                       | S/P/WB                 |
| SAA                           | Bacterial/Virus Infection                         | S/P/WB/Fingertip Blood |
| SAA/CRP                       | Neonatal Sepsis, Bacterial/Virus Infection        | S/P/WB/Capillary Blood |
| <b>Diabetes Mellitus</b>      |                                                   |                        |
| HbA1c                         | Diabetes Mellitus                                 | WB                     |
| <b>Reproduction/Fertility</b> |                                                   |                        |
| FSH                           | PCOS, Infertility Evaluation, Pituitary Disorders | S/P/WB                 |
| LH                            | PCOS, Infertility Evaluation                      | S/P/WB                 |
| Prog/HCG                      | Pregnancy Test                                    | S/P/WB                 |
| <b>Thyroid Function</b>       |                                                   |                        |
| T3                            | Hyperthyroidism, Hypothyroidism                   | S/P/WB                 |
| T4                            | Hyperthyroidism, Hypothyroidism                   | S/P/WB                 |
| <b>Metabolic Marker</b>       |                                                   |                        |
| 25-OH-VD                      | Osteomalacia, osteoporosis                        | S/P/WB                 |

| Test Items                             | Dry Immunofluorescence | Dry Clinical Chemistry  | Sample Types           |
|----------------------------------------|------------------------|-------------------------|------------------------|
| GLU/HbA1c                              | HbA1c                  | GLU                     | WB/ Fingertip Blood    |
| Uric Acid                              |                        | Uric Acid               | S/P/WB/Fingertip Blood |
| TCH                                    |                        | TCH                     | S/P/WB/Fingertip Blood |
| TG                                     |                        | TG                      | S/P/WB/Fingertip Blood |
| HDL-C                                  |                        | HDL-C                   | S/P/WB/Fingertip Blood |
| Lipid profile                          |                        | Lipid profile           | S/P/WB/Fingertip Blood |
| Lipid profile/Uric acid                |                        | Lipid profile/Uric acid | S/P/WB/Fingertip Blood |
| Lipid profile/Uric acid/CK-MB/cTnI/Myo | CK-MB/cTnI/Myo         | Lipid profile/Uric acid | S/P/WB/Fingertip Blood |
| Cr                                     |                        | Cr                      | S/P/WB/Fingertip Blood |
| Cr/mAlb                                | mAlb                   | Cr                      | Urine                  |
| Urea                                   |                        | Urea                    | S/P/WB/Fingertip Blood |
| Urea/Cr                                |                        | Urea/Cr                 | S/P/WB/Fingertip Blood |
| ALT                                    |                        | ALT                     | S/P/WB/Fingertip Blood |
| AST                                    |                        | AST                     | S/P/WB/Fingertip Blood |
| ALT/AST                                |                        | ALT/AST                 | S/P/WB/Fingertip Blood |
| ALT/AST/Urea/Cr                        |                        | ALT/AST/Urea/Cr         | S/P/WB/Fingertip Blood |
| HCY/CK-MB/cTnI/Myo                     | CK-MB/cTnI/Myo         | HCY                     | S/P/WB/Fingertip Blood |



Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China  
 Tel: +86-25-68568508/68568594  
 Fax: +86-25-68568500  
 E-mail: sales@getein.com.cn; overseas@getein.com.cn  
 Web: www.getein.com



FSC CE NMPA IVD

