

## AP-100

Automated Dosing System

An automated and compact design for accurate and safe dosing of Mindray ID/AST




# AP-100

## Automated Dosing System



The multi-location technology helps to reduce manual intervention and improve both accuracy and speed.



The compact design allows for plate dosing in the biosafety cabinet, which ensures higher safety during operation.

### Technical specifications

Sampling principle	The syringe and pump work together for aspirating and dispensing
Display	LED digital screen for displaying test card categories
Sampling range	TDR YEAST-96, TDR STAPH-96, TDR STR-96, TDR NF-96, TDR ONE-96, TDR NH-96
Throughput	30 test cards per hour
Reproducibility	CV ≤ 5%
Power supply	AC100~240V, 50/60H
Dimension and weight	340mm×300mm×360mm 13±1kg
Working conditions	Ambient temperature: 15~30℃, Relative humidity: 20%RH~80%RH, no condensation Altitude: -400~3000m

[www.mindray.com](http://www.mindray.com)

P/N:ENG-AP-100-210285X2P-20221020

©2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

**mindray**  
healthcare within reach